

Certified ISO 13485:2016, CE cGMP / FDA Approved



Concise[®]
Enterprises

(Surgical Instruments Manufacturing Company)



DIAGNOSTIC INSTRUMENTS CATALOGUE

- DIAGNOSTIC E.N.T
- LARYNGOSCOPE
- OTOSCOPE
- TONGUE DEPRESSERS,
- HEAD LIGHTS, HUMMERS

Concise Enterprise the name of ergonomics, long-lasting quality and steady progress, was established in 1958 as a Manufacturers **Surgical Instruments** **Dental Instruments, Scissors.**

Certification :-

Concise Enterprise is certified with Quality management system ISO 9001:2008, ISO 13485:2016 CE Mark, cGMP (F.D.A / USA) approved.

Health and Safety Policy, Environment policy always adopted as per international standard.

Quality policy:-

- Concise is Working for Customer satisfaction
- Customer Confidence
- Produce cost effective quality products
- Doing things right at first time.
- Like the partner that is behind you.

Founder Mr. Mohammad Yousaf (Late) devoted his life to serve the world in surgical / dental field. Following the footsteps of their fore-fathers, they have engaged themselves in the manufacturing instruments made of Stainless Steel. The experience and knowledge transferred from the fore-fathers to founder count to be the back bone of excellent quality **CONCISE** products.

Organizational Structure :- We are trying our best to provide you the best service and atmosphere to do business. Today, **Concise Enterprises** well known throughout world for its workmanship, Prime Quality, Competitive Prices, and Prompt Services. **Concise** is equipped with modern facilities to manufacture Instruments according to needs of its valued customer.

Concise Enterprises is thankful to all of its worthy customers for their trust and for their appreciation of quality and services **Concise** provide them.

Regards,
Concise Team

Certificates





DIAGNOSTIC SET CED-41-01001

2.5 V Standard illumination for use with
2 batteries C type.

Including:

Otoscope Head: Luir Lock fitting.
Chrome plated.

Fix lens 4 x in plastic mount can be
removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luir Lock fitting.

Lens selection wheel permits choice of corrective
lenses from 0 to - 20 and 0 to + 20 diopters.

Single large circle aperture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in our new beautiful
look and quality Blue Case, which protect our quality
Products from moisturizing and damaging.



With 2 spare Bulbs





DIAGNOSTIC SET

CED-41-01002

2.5 V Standard illumination for use with 2 batteries type C.

Including:

Otoscope Head: Luir Lock fitting.

Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: pin contact fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to + 20 diopters.

Single large circle operture.

Bent-arm illuminator: Luir Lock fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

The Diagnostic set is supplied in new beautiful look and quality Blue Case, which protect our quality

Products from moisturizing and damaging.





DIAGNOSTIC SET

CED-41-01003

2.5 V Standard illumination for use with 2 batteries C type.

Including:

Otoscope Head: Luir Lock fitting. Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luir Lock fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to + 20 diopters.

Single large circle operture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in new beautiful look and quality Red Wooden Case, which protect our quality Products from moisturizing and damaging.



DIAGNOSTIC SET

CED-41-01004

2.5 V Standard illumination for use with 2 batteries C type.

Including:

Otoscope Head: Luir Lock fitting. Chrome plated.

Fix lens 4 x in plastic mount can be removed for introduction of instruments.

Specula sizes: 2.5, 3.5 & 4.5 mm Included,

Ophthalmoscope head: Luir Lock fitting.

Lens selection wheel permits choice of corrective lenses from 0 to - 20 and 0 to + 20 diopters.

Single large circle operture.

Bent-arm illuminator: pin contact fitting.

Nasal speculum: Chrome plated.

Two laryngeal mirrors: n. 4 and 5.

Battery handle: 28mm. For 2 cell C-batteries.

Chrome plated. For adjustment of all illumination.

Replacement bulb: 2.5 V

Diagnostic set is supplied in new beautiful look and quality Black Case, which protect our quality Products from moisturizing and damaging.



MOUTH MIRRORS DIAGNOSTIC SET

CED-41-01005

10 Disposable Mirrors in plastic fitting.

01 Mirror Handle with Led Bulb having inside 01 set of 3 button cells each button cell hare 1.5 V.

This set here complete 03 sets of button cells

Packed in beautiful Transparent Plastic case.





DIAGNOSTIC PENTOSCOPE SET CED-41-01006

2.5 V Standard illumination for use with 2 batteries type AA.

Consisting:

Universal Handle thread type fit for all Pentoscope Heads like.

- Otoscope Head.
- Ophthalmoscope Head.
- Throat Branch Head.
- Mouth Mirrors Head.

It have two exploring mirrors No. 4 & 5 mm.

Its spare 06 disposable Plastic ear Canulas are with sizes 2 x 2.50 mm, 2 x 3.50 mm & 2 x 4.5 mm.

On special order it can be as 2 x 1.5 mm for new born babies 2 x 2.5 mm for pediatric and 2 x 4 mm for adults, so it means total 06 canula also available in a plastic tube packing.

2 spare 2.5 V bulb.


Pentoscope set battery handle is operating with two AA size batteries.

This Pentoscope set Packed in our new beautiful Blue Colour plastic case in same blue inner packing material.

Spare parts also available of any item.



MC-INTOSH FIBER OPTIC LARYNGOSCOPE BLADE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Mc-Intosh Laryngoscope Fiber Optic

Halogen illumination. Manufactured from AISI 304 Stainless-Steel. Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Laryngoscope Battery
Handle for Adult C type 2 batteries

Available Sets with Battery Handle

CED-42-02001

Set of 0,1,2 (**Pediatric set**)
Supplied in Green plastic case

CED-42-02002

Set of 1,2,3 (**Standard set**)
Supplied in Green plastic case


CED-42-02003

Set of 1,2,3,4 (**Adult set**)
Supplied in Green plastic case



Mc-Intosh Fiber Optic Laryngoscope Set
Halogen Illumination
Manufactured by AISI: 304 Stainless Steel Material

MC-INTOSH LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Fiber Optic Laryngoscope

Plastic case, Halogen illumination. Manufactured from AISI 304 Stainless Steel Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

Blades Availability:

- * Matt Finish
- * Mirror Finish
- * Sand Blast Finish
- * Pls, confirm in your order list

Available Sets with Battery Handle

CED-42-02004

Set of 0,1,2 (Pediatric Set)
Supplied in Green plastic case


CED-42-02005

Set of 1,2,3 (Standard Set)
Supplied in Green plastic case

CED-42-02006

Set of 1,2,3,4 (Adult Set)
Supplied in Green plastic case

Fiber Optic Laryngoscope Miller



BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

Available Sets with Battery Handle

CED-42-02007

Set of 0,1,2 (Pediatric Set)
Supplied in Green plastic case

CED-42-02008

Set of 1,2,3 (Standard Set)
Supplied in Green plastic case

CED-42-02009

Set of 1,2,3,4 (Adult Set)
Supplied in Green plastic case



Laryngoscope Battery Handle for Adult C type 2 Batteries



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Mc-Intosh Laryngoscope Set of Fiber Optic Halogen Illumination
Manufactured From: AISI 304 Stainless Steel Material



Laryngoscope Battery Handle for Adult C type 2 Batteries



Fig. 1



Fig. 2




Fig. 3



Fig. 4

MC-INTOSH LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

Available Sets with Battery Handle

CED-42-02010

Set of 0,1,2 (Pediatric Set)
Supplied in Blue plastic case

CED-42-02011

Set of 1,2,3 (Standard Set)
Supplied in Blue plastic case

CED-42-02012

Set of 1,2,3,4 (Adult Set)
Supplied in Blue plastic case

Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

Blades Availability:

- * Matt Finish
- * Mirror Finish
- * Sand Blast Finish
- * Pls, confirm in your order list



With two Spare Bulbs



Laryngoscope Battery
Handle for Adult C type 2 batteries



Fig. 1



Fig. 2




Fig. 3



Mc-Intosh Laryngoscope Conventional
Standard Illumination
Manufactured From Stainless Steel material

MC-INTOSH LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

Available Sets with Battery Handle

CED-42-02013

Set of 0,1,2 (Pediatric Set)
Supplied in Gray plastic case

CED-42-02014

Set of 1,2,3 (Standard Set)
Supplied in Gray plastic case

CED-42-02015

Set of 1,2,3,4 (Adult Set)
Supplied in Gray plastic case

Mc-Intosh Laryngoscope


Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel Matt Type, Standard Illumination 2.5 V. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Gray plastic case.

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



Laryngoscope Miller



BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

Available Sets with Battery Handle

CED-42-02016

Set of 0,1,2 (Pediatric Set)
Supplied in Blue plastic case

CED-42-02017

Set of 1,2,3 (Standard Set)
Supplied in Blue plastic case

CED-42-02018

Set of 1,2,3,4 (Adult Set)
Supplied in Blue plastic case

Mc-Intosh Laryngoscope

Plastic case, Standard illumination. Manufactured from AISI 304 Stainless Steel. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Blue plastic case.

Blades Availability:

- * Matt Finish
- * Mirror Finish
- * Sand Blast Finish
- * Pls, confirm in your order list



With two Spare Bulbs



Laryngoscope Battery Handle for Adult C type
2 batteries



Fig. 1



Fig. 2



Fig. 3



Fig. 4



Laryngoscope Set of conventional Miller
Standard Illumination
Manufactured From Stainless Steel Material

Conventional Set



Available blade sizes 2,3 & 4 with
'AA' & 'C' type battery handles

'C' type handle

CED-42-02019

Flexible Laryngoscope Conventional

Halogen Illumination

Manufactured From Stainless Steel High grade
304 material.

Fiber Optic Set



CED-42-02020

Flexible Laryngoscope Set of Fiber Optic

Halogen Illumination

Manufactured From Stainless Steel High grade
304 material packed in green plastic box.




Available Blade sizes 2,3 & 4 with
'AA' & 'C' type battery handles

'AA' type handle

Laryngoscope Battery
Handle for Adult AA Batteries

Conventional Set



Blades

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm




CED-42-02021
Conventional Set,
Mac + Miller Combo Double Set
of 10 Pcs.

4 Miller Blades
4 Mac Blades
1 Handle Small
1 Handle Medium

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list

MC-INTOSH LARYNGOSCOPE



Blades


Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm

Mc-Intosh Laryngoscope Conventional

Standard illumination.

Manufactured from AISI 304 Stainless Steel
All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Wooden made case.

Fiber Optic Set



Blades

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm




CED-42-02022
Fiber Optic
Mac + Miller Combo Double Set
of 10 Pcs.

4 Miller Blades
4 Mac Blades
1 Handle Small
1 Handle Medium

Blades Availability:

- ★ Matt Finish
- ★ Mirror Finish
- ★ Sand Blast Finish
- ★ Pls, confirm in your order list



Blades

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Mc-Intosh Laryngoscope
Fiber Optic

Halogen illumination. Manufactured from AISI 304 Stainless Steel Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green Wooden Made case.



**Standard Otoscope
CED-43-03001**

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luer / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Standard Otoscope
CED-43-03002**

This otoscope is 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle is Pin Serew lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Case Packing



**Standard Otoscope
CED-43-03003**

This otoscope is 2.5 V standard illumination, handle is powder coated as well otoscope head too. Handle and head is Pin Screw lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Nice Wooden Box Packing



**Mini Otoscope Fiber Optic
CED-43-03004**

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Lonized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colours.



**Mini Otoscope Fiber Optic
CED-43-03005**

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Lonized Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v the diagnostic set is supplied in Black Pouch. Available in multi colours.



**Standard Otoscope + Ophthalmoscope Set
CED-43-03006**

This Set is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Plastic Box Packing.



'C' type
Battery handle



**Standard Ophthalmoscope Set
CED-43-03007**

A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image Packed in Black Plastic Case.



'C' type Battery Handle



**Standard Otoscope Set
CED-43-03008**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing



**Standard Otoscope Set
CED-43-03009**

this otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing




**Otoscope + Ophthalmoscope Set
CED-43-03010**

This otoscope is with 2.5 V standard illumination, handle is chrome plated as well otoscope head too. Handle and head is Luir / Faster lock fitting. Head of otoscope is with 6x magnifying glass. 03 Resuable speculas with sizes 2.5, 3.5 & 4.5 mm. A Full Functional Ophthalmoscope compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image A Full Functional Ophthalmoscope Battery handle for 2 'c' batteries. Rheostat for adjustment of illumination. Supply in Black Pouch Packing.



**Mini Ophthalmoscope
CED-43-03011**

A Full Functional Ophthalmoscope.
compact durable, light weight.
Standard bulb and precision system
4 different apertures selection
range of lenses -20 to +20
bright white light and brilliant image
Made of plastic Material.



**Mini Ophthalmoscope
CED-43-03012**

A Full Functional Ophthalmoscope.
compact durable, light weight.
Standard bulb and precision system
4 different apertures selection
range of lenses -20 to +20
bright white light and brilliant image
made of plastic Material.



**Standard Ophthalmoscope
CED-43-03013**

A Full Functional Ophthalmoscope.
compact durable, Standard bulb and
precision system
4 different apertures selection
range of lenses -20 to +20
bright white light and brilliant image
Handle made of Stainless Steel Material.



Mini Ophthalmoscope


CED-43-03014 - CED-43-03019

A Full Functional Ophthalmoscope. compact durable, light weight. Standard bulb and precision system 4 different apertures selection range of lenses -20 to +20 bright white light and brilliant image. Available in multi colours.

Otoscope Fiber Optic

CED-43-03020 - CED-43-03025

2.5 V fiber optic halogen light. Otoscope head: pin contact fitting. Fixed lens 6X in plastic mount may be removed for introduction of instruments. Reusable specula 2, 3 and 4 mm. Included Battery Handle Made of Brass, Lonzed Plated handle with Push button for 2 AA Batteries. Replacement bulb: 2.5 v Available in multi colours.



Otoscope Fiber

CED-43-03014
CED-43-03015
CED-43-03016
CED-43-03017
CED-43-03018
CED-43-03019

Ophthalmoscope

CED-43-03020
CED-43-03021
CED-43-03022
CED-43-03023
CED-43-03024
CED-43-03025



(A Set of Four Blades)

**Tongue Depressor Holder
CED-44-04001**

2,5 V standard illumination. Two adult, two child tongue blade holders. 304 stainless steel blade easy to auto clave 2 AA plastic handle in six different colours. A to F

**Tongue Depressor Holder
CED-44-04002**


2,5 V standard illumination. Two adult, two child tongue blade holders. 304 stainless steel blade easy to auto clave 2 AA chrome plated handle with push ON/OFF system.



**Tongue Depressor
Holder Plastic Body
CED-44-04003**



**Tongue Depressor
Holder metal Body
CED-44-04004**



**Tongue Depressor
Conventional**

**Tongue Depressor Holder
CED-44-04005**

2,5 V standard illumination. Two adult, two child tongue blade holders. 304 stainless steel blade easy to auto clave 2 AA chrome plated handle with push ON/OFF system.



Head Light
CED-44-04006

Head light with standard illumination 2,5 V.
Chrome plated head with magnifying glass.
Movable in all directions by means of a joint.
Head band can be adjusted by metal screw.
Supplied with cable for direct connections to the battery handle.

Battery Handle: Chrome plated for 2 C batteries.
Battery Handle: Plastic made for 2 AA batteries.



Head Light
CED-44-04007

Deluxe head light with standard illumination 2,5 V.
Plated head with magnifying glass.
Movable in all directions by means of a joint.
Padded Head band can be adjusted by metal screw.
Supplied with cable for direct connections to the battery.

Battery Handle: Chrome plated for 2 C batteries.
Battery Handle: Plastic made for 2 AA batteries.



Head Mirror
CED-44-04008

The black plastic head band can be adjusted to the individual head size by means of the metal buckle.
The hinge allows adjustment of the mirror, Ø mm. 90, in all direction.



Light For Speculum
CED-44-04009

Light for Cusco and Graeve speculum.
Battery Handle: Made in plastic material.
19mm, Ø for 2 AA batteries,
Cable with bulb and clip for fixation.
Standard illumination 2,5 V.

Battery Handle: Chrome plated for 2 C batteries.
Battery Handle: Plastic made for 2 AA batteries.



CED-44-04010

Varioflex combination hammer according to Fassbender, with never needle and brush. Solid handle with centimeter scale marking (approx. 155g) length 20 cm.

CED-44-04011

Tromner hammer original model heavy version (approx. 210g)



CED-44-04012

Rossier - Reflex hammer, with elastic handle, length: 36cm (approx. 95g)



CED-44-04013

Rossier - Reflex hammer, with elastic handle, length: 24cm (approx. 90g)



CED-44-04014

Babinski - Reflex and percussion hammer, interchangeable with screw fitting, with metal grip (approx. 95g)



CED-44-04015

Babinski - Reflex and percussion hammer, interchangeable with screw fitting, with metal grip (approx. 90g)



CED-44-04016

Buck, Neurological Reflex hammer, with metal handle, brush and needle (approx. 90g)



A



B



C




D

CED-44-04017

Colorflex, reflex and percussion hammer, Small version, length: 17cm
Color: red, blue, yellow, black (approx. 100g)

Conventional

MC-INTOSH LARYNGOSCOPE



Blades

Size	A	B	C	Working Length
0	80mm	58mm	30mm	55mm
1	97mm	76mm	36mm	70mm
2	111mm	88mm	36mm	78mm
3	133mm	112mm	36mm	108mm
4	154mm	131mm	38mm	133mm



'AA' Battery
Handle



Fig. 1



Fig. 2



Fig. 3




Fig. 4

Mc-Intosh Laryngoscope Thread Type

Rexien Pouch, Standard illumination. Manufactured Chrome plated handle and S.S Made Blades Thread type. Standard Illumination 2.5 V. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in plastic Pouch.

Laryngoscope Miller
Thread Type in Rexien Pouch



Blades

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm




CED-45-05001
Mac Economy Set
Set of 4 Blade + 1 Handle



CED-45-05002
Miller Economy Set
Set of 4 Blade + 1 Handle

MC-INTOSH FIBER OPTIC LARYNGOSCOPE



BLADE

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

Mc-Intosh Laryngoscope Fiber Optic Set


Halogen illumination. Manufactured from AISI 304 Stainless-Steel, Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic economy case.

Blades Availability:

- * Matt Finish
- * Mirror Finish
- * Sand Blast Finish
- * Pls, confirm in your order list

Fiber Optic Miller Laryngoscope

Thread Type in Rexien Pouch



Blades

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm

Fig. 1

Fig. 2

Fig. 3

CED-45-05003

Mc-Intosh Economy Set
Set of 3 Blade + 1 Handle

Fig. 1

Fig. 2


Fig. 3

Fig. 4

CED-45-05004

Miller Economy Set
Set of 4 Blade + 1 Handle

Bulbs For Conrentional + Fiber Optic Use



CED-45-05005 Standard Bulbs	CED-45-05006 Frosted CED-45-05007 Clear 1-8 72NS3 Thread 2.5V .28A 20 Hour Life Nickel Plated Vacumm T-1 Lamp	CED-45-05008 Frosted CED-45-05009 Clear 8-32 UNC 2A Thread 2.5V .28A 20 Hour Life Nickel Plated Vacumm TL-1 ½ Lamp	CED-45-05010 8-32 UNC 2A Thread 2.5V .38A 20 Hour Life Nickel Plated Vacumm TL-1 ½ Lamp	CED-45-05011 8-32 UNC 2A Thread 2.5V .35A 15 Hour Life Nickel Plated Vacumm TL-1 ¼ Lamp	CED-45-05012 3.5 X .35mm Thread 2.5V .28A 15 Hour Life Nickel Plated Vacumm TL-1 ¼ Lamp	CED-45-05013 8-32 UNC 2A Thread 2.5V .30A 15 Hour Life Nickel Plated Vacumm TL-1 ¼ Lamp
--	--	---	--	--	--	--



Fiber Optic Handles With Bulb

	Size	Batteries
CED-45-05014	Medium / Standard	2C
CED-45-05015	Small / Penlight	2AA
CED-45-05016	Stubby	2AA
CED-45-05017	Large	2D



Conventional Handles

	Size	Batteries
CED-45-05018	Medium / Standard	2C
CED-45-05019	Small / Penlight	2AA
CED-45-05020	Stubby	2AA
CED-45-05021	Large	2D



CED-45-05022
Disposable specula LTI-1
Size: 1.5, 2.5, 3.5, 4 & 4.5mm




CED-45-05023
Reusable specula LTI-2
Size: 3.5, 4.5 & 5.5mm



CED-45-05024
Reusable specula LTI-3
Size: 2.5, 3.5 & 4.5mm

FIBER OPTIC MILLER LARYNGOSCOPE



BLADE

Size	A	B	C
0	80mm	55mm	38mm
1	97mm	75mm	40mm
2	155mm	131mm	40mm
3	194mm	170mm	40mm
4	256mm	231mm	48mm
5	300mm	271mm	48mm

Available Sets with Battery Handle

CED-45-05025

Set of 0,1,2
Supplied in Black Wooden case

CED-45-05026

Set of 1,2,3
Supplied in Black Wooden case

CED-45-05027

Set of 0,1,2,3,4,5
Supplied in Black Wooden case

Fig. 0

Fig. 1

Fig. 2

Fig. 3

Fig. 4

Fig. 5



Otoscope

2.5 V Standard light.
Otoscope head: Luer Lock contact fitting.
Fixed lens 4X in plastic mount may be removed for introduction of instruments.
Reusable Longer veterinary speculas 2.5, 3.5 & 4.5mm dia. Included Battery Handle Made of Brass, Ionized Plated handle with Push button for 2 'C' Batteries.
Replacement bulb: 2.5 v



Laryngoscope Battery for 2 'C' type batteries



Miller Laryngoscope Fiber Optic CED-45-05028

Halogen illumination. Blades manufactured from stainless steel AISI 304.
Supplied in Red / Black wooden Box.



*Unique In
Quality*



Concise
Enterprises

(Surgical Instruments Manufacturing Company)

Ahmed Din Estate, Noul More, Roras Road,
Sialkot-51310, Pakistan

+92 300 9616641 | Tel. +92 52 3574077 | Fax +92 52 3574078

E-mail: inquiry@concise.com.pk
export@concise.com.pk
contact@concise.com.pk
infoconcise@gmail.com

<http://www.concise.com.pk> | www.medical-instruments.biz

USA Office: Concise Med & Beauty USA

Astoria NY 11105-5746

Tel: +1 (347) 3050081

E-Mail: usa@medical-instruments.biz



Certificates

